

UV Light Box - 110V

MYSTAIRE®



Made in the USA



GENERAL SPECIFICATIONS

Construction Material	Work Surface	Polypropylene
	Chamber	Polypropylene
	Front Sash	Clear Polycarbonate
	Head Assembly	Polypropylene
Controls		Micro-processor Controller
UV Lighting		254-nm Shortwave UV
UV Monitoring Package		Audible / Visible Alarm
UV Safety Switch		Included
Automatic UV Timer		User Defined Timed Irradiation
Lab Event Timer		Adjustable Audible / Visible Timer

UNIT DIMENSIONS*

Product #	External Dimensions			Internal Dimensions			Door Dimensions	
	Width	Depth	Height	Width	Depth	Height	Width	Height
MY-DB24C	24"	24"	15"	22"	23"	7"	23"	9"
MY-DB42C	42"	24"	15"	40"	23"	7"	41"	9"

*All dimensions are nominal within +/- 0.5"

TECHNICAL SPECIFICATIONS

Product #	Electrical		Unit Weight**
	Voltage	Power	
MY-DB24C	110 V	70 W	150 lbs.
MY-DB42C	110 V	90 W	195 lbs.

**All weights are an approximation

www.mystaire.com

Phone: +1 (919) 229-8511 • Toll Free: (877) 328-3912 • Email: contact@mystaire.com